

Professional Hearing Solutions

Latest Hearing Aids



DATA SHEET

www.professionalhearingsolution.com 0332-5014111




BTE HEARING SYSTEMS HP G4

Tech Level 3



BATTERY: 675 AMPLIFICATION: 84 dB



HP G4



BASIC FEATURES

	Tech Level	3
Earhook		•
Battery compartment as on/off switch		•
Programmable push button		•
Programmable rocker switch		•
Level-dependent signal tones/melodies (can be activated/deactivated)		•
Can be used with either disposable or recharge	able battery	_
Telephone coil		•
Status LED, programmable		_
Personal color concept		-
Audiomatic power-on delay (can be activated/deactivated)		•
IP67-certified		•

SIGNAL PROCESSING

Anti-Feedback system	•
Noise Manager	
> Adaptive noise reduction	😑 (on/off)
> Wiener filter	😑 (on/off)
> Adaptive wind noise reduction, binaural	-
> Adaptive wind noise reduction	-
> Impulse suppressor	_
> Auto-situation adaption	-
> Selectronic	_
AudioTronic multi-microphone system	
> Panorama	-
> Directional static	•
> Automatic	•
> Adaptive	_
> Speech 360	-
> AudioFocus 360	_
> AudioSpot	_
Frequency and dynamics concept	
> TRC S	•
> Selective frequency compression	-
> HiFi functionality	_
> Sound Upgrade	_
(in Wireless Streaming mode)	
> Gain setting dependent on environment	_
(only in Wireless Audio Streaming mode)	-
Programmable tinnitus function	_
AUTOMATIC FUNCTIONS	

Comfort365-Acclimatic-Intelligent Acclimatic-Comformatic-

TECHNICAL FEATURES

	Tech Level	3
Signal processing channels		12
Frequency channels		6
AGC channels		6
MPO channels		6
Hearing programs		4
> MusicSelect		-
> 2earPhone		-
> EchoClear/dereverberation		-
Data Logging		•
Wireless		
> AudioLink		_
> Binaural synchronization		-
> Wireless audio streaming		-
> CROS/BiCROS		-

ACCESSORIES | OPTIONS

Clip for battery compartment in red and blue for side recognition		•
Thin Tube with Open Tip		_
Smart Power		-
Audio streamer Smart Connect	8	_
Smart Remote remote control		-
Remote control Smart Key	8	0
Small earhook		0
AudioFix		_
Fitting Set - Thin Tube		-
Audio shoe set		0
Battery compartment safety		_
AutoPhone set		-
APPS		
APPS Smart Connect app	6	_
	ୖ	
Smart Connect app	8	 0
Smart Connect app		- - 0 -
Smart Connect app with AudioDirSelect Smart Remote app		- - 0 -
Smart Connect app with AudioDirSelect Smart Remote app		- - 0 -
Smart Connect app > with AudioDirSelect Smart Remote app > with AudioDirSelect		- - 0 -
Smart Connect app > with AudioDirSelect Smart Remote app > with AudioDirSelect PROGRAMMING		- - 0 -
Smart Connect app > with AudioDirSelect Smart Remote app > with AudioDirSelect PROGRAMMING Programming adapter 312		- - - -


• = Standard equipment O = optional — = not available



HP G4 with earhook


MAXIMUM POWER OUTPUT

LE = 90 dB



MAXIMUM AMPLIFICATION

Amplification at LE = 50 dB



dB 70 100 Hz 1000 Hz 10000 Hz 1000 Hz 1000 Hz 10000 Hz

TECHNICAL INFORMATION

140 dB	144 dB
130 dB	137 dB
132 dB	138 dB
80 dB	84 dB
65 dB	71 dB
69 dB	73 dB
56 dB	64 dB
675	675
184	184
100 – 5,600 Hz	100 – 5,900 Hz
2.4 mA	2.4 mA
24 dB	26 dB
_	
102	107
3% 2% 1%	7% 3% 2%
	130 dB 132 dB 80 dB 65 dB 69 dB 56 dB 675 184 100 – 5,600 Hz 2.4 mA 24 dB — 102 3% 2%

 $^{\eta}$ Technical data measured in accordance with IEC 60118-7:2005 and ANSI S3.22-2009 at 2 ccm coupler $^{2)}$ Technical data measured in accordance with IEC 60118-0 at ear simulator

 ▲ WARNING Small parts present a choking hazard. This device is not suitable for fitting to babies, small children or mentally disabled persons.
 ▲ WARNING The maximum output sound pressure level of the hearing systems can reach or exceed 132 dB SPL. Risk of damage to the hearing of the wearer. Ensure that the hearing systems are fitted with care.

For control ranges and more programming features see Hearing System Simulation of Connexx 8.4, AudioFit 8.4.5 or higher.

tion